

# NMR-Al powered Drug Discovery



CEEC Mee 2025

Prague 09.09.2025

www.aiffinity.com

info@aiffinity.com

Purkyňova 127, 612 00 Brno-Medlánky, Czech Republic

## Our solution: NMR-Al Integrated platform for Streamlined Target-to-Lead Process



Examining how drugs bind to the target protein



Finding drug candidates "hits" through experiments & computations



Refining best hits to become lead drug candidates













3D & 4D NMR

1D NMR

1D & 2D NMR







## We Scrutinize Atomic Interactions with NMR

- We measure both sides:
  - $\circ$  2D NMR (protein)  $\rightarrow$  dozens of data points/day
  - o ultra-fast 1D NMR (ligand) → hundreds of data points/day
- Our NMR-Al platform fuses these data to discover new hits, grow/link fragments, and iterate design









## deepHitExplorer

Multi-objective Neural Network deepHitExplorer Virtual Input Genetic Algorithm Screening Viewer Molecule libraries 1D 1H Reference spectrum Not-binding 41 40 39 38 37 36 35 34 33 32 31 30 29 Ligand Epitope Compounds Al tool integrated with NMR Diverse molecules Ligand-observed Interactive exploration



NMR spectrum

of the diverse

chemical space

data for enhanced molecule

scoring

search based on

chemical structure,

predicted score and solubility

undergoing

epitope prediction

and scoring

## deepHitExplorer Viewer

Fix selection and let the program suggest other chemically diverse, top-scored, soluble molecules

Filtering based on the molecule properties


Dive in the 3D scatter plot selecting the desired molecules for experimental testing

Compare the picked molecules



Pick and download your selection

Visualize
the
predicted
ligand
epitope
and
molecular
properties





www.aiffinity.com

#### **Protein NMR**

#### **NMR**

#### **NMR**

#### X-ray Cryo-EM **NMR AlphaFold**

#### Spectral


Resolution: 2D < 3D < 4D < 5D

Sensitivity: 2D > 3D > 4D > 5D

**Total time:** 2D < 3D < 4D < 5D

**Complexity:** 2D < 3D < 4D < 5D







protein









#### Our competitive advantages:

- Large proteins (>20 kDa)
- Intrinsically Disordered Proteins






### Interaction mapping by NMR – Backbone assignment



Protein backbone resonance assignment

5 x 3D NMR spectra




Backbone & sidechain resonance assignment + 3D structure

5 x 3D NMR spectra + 2 x 4D spectra



## 4D-GraFID: Superior Backbone Assignment






www.aiffinity.com

## Androgen Receptor (AR-V7)

- IDPs (>30% of human proteome) lack defined 3D structure
- Ligand binding reweights the ensemble unpredictably → structurebased approaches are unreliable
- What matters: atom-level contact maps (protein atoms ↔ ligand atoms)
- Our drug design tech in based on atomic interactions captured by NMR both from the ligand's and the protein's perspective → ideal for IDPs



conformational ensemble from *Metadynamic Metainference J. Chem. Theory Comput. 2025, 21, 9, 4898–4909* 



### 4D-GraFID: protein structure from fewer spectra


Chemical shift assignment (**4D-GraFID**)



4D HCNH NOESY 4D HCCH NOESY

#### **Backbone spectra:**

3D CBCA(CO)NH 3D HNCACB 3D HN(CO)CA 3D HNCA 3D HNCO




NMR restraint optimization

3D structure model generation

powered by






Final model

Medium-range and long-range NOE distance restraints






- **Modified OpenFold**
- nElt (248 amino acids)
- **Only 4D HCNH NOESY** spectrum


1.8 mM, 20 mM sodium phosphate (pH 6.5), 100 mM NaCl, 5% D2O, 37°C.


Video on **YoutTube** 




## 4D-GraFID: protein structure from 2x 4D spectra







## SaaS platform available online or on-premise





## Our Wet-Lab Facilities for Protein Production & Biochemical Assays




A) General bench & centrifuge

B) Class II biosafety cabinet




C) Plate reader & liquid handling



D) Automation & PCR)

#### Privileged access to <u>full-fledged biophysical facilities</u>










- structure determination (NMR, X-ray, SAXS)
- compound screening (MST, SPR)
- binding validation & interactions mapping (HDX-MS, CD, SEC, ITC)


## Joint Drug Discovery powered by our NMR-Al tech





#### R&D Team

#### Administrative Team





Research & Development Founder, CEO & CTO



M.Eng. Research & Development Computational Biologist



Administration Backoffice & Grant Manager



Ekaterina Burakova, M.Sc., Dr.rer.nat.



Research & Development NMR & Data Scientist



Research & Development Computational Chemist

SaaS Dev Team

Yevgeniy Fomin



Software Developer



Valerij Talagayev, M.Sc.



Research & Development Chemoinformatics Developer



Research & Development Protein Biochemist

**Affiliated Contributors** 





Research & Development **Application Scientist** & Project Manager



Research & Development PhD Student

Václav Hanzl, Ing., CSc.



ML Expert

#### Manolis Arsenikos, M.Pharm.



Research & Development PhD Student



Research & Development Protein Biochemist

Tereza Vučková, M.Sc., Ph.D.



Research & Development Protein Biochemist



#### "Focus and excel"

#### Al | ffinity s.r.o.

IČO 11992719 Purkyňova 127, 612 00 Brno-Medlánky, Czech Republic

www.aiffinity.com

info@aiffinity.com

## **Open Positions:**

#### **Cheminformatics Developer**

https://www.linkedin.com/jobs/view/4293495038/

#### PhD position in Cheminformatics-NMR

https://www.linkedin.com/jobs/view/4296103102/

#### **Investors:**







### 1D NMR methods in Drug Design

N≡XMR

- Photo-CIDNP, STD, WaterLOGSY
- Photo-CIDNP: more sensitive, faster (1,500 samples/day), benchtop NMR
- Protein 3D structure not necessary
- Give the atoms of the small molecule that interact with the protein (ligand epitope)
- We develop both Generative AI and Virtual Screening solutions incorporating ligand epitope information

